

Vassar Diagnostics Supply Order Form

Practice Name: _____

Contact Name: _____

Contact Phone Number: _____

Specimen Collection		
QTY	UNITS	DESCRIPTION
	100/Flat	Blue Sodium Citrate Tubes
	100/Flat	Gray Sodium Fluoride Tubes
	100/Flat	Green Lithium Heparin Tubes - 4m1
	100/Flat	Lavender EDTA Tubes - 3m1
	100/Flat	Pink EDTA Tubes -10m1
	100/Flat	Red top tube - 3m1
	Flat	Gold Top Tubes - 5m1
	Flat	Gold Top Tubes - 3.5m1
	Box	Vacutainer 22G (Black) Vacutainer needle
	Box	Vacutainer 21G (Green) Vacutainer needle
	Each	Vacutainer plastic holder (hub)
Miscellaneous		
QTY	UNITS	DESCRIPTION
	Each	24hr urine collection container
	Each	Sterile Container
	Each	Plastic pathology container
	Pack	Specimen Transport Bag
	Pack	STAT Specimen Transport Bag
	Each	Specimen Transport Blue Frozen Bag
	Each	Specimen Box
	Each	Centrifuge
	Package	Tourniquet 1X18" Latex Free
Microbiology		
QTY	UNITS	DESCRIPTION
	Each	Anerobic Blood Culture Bottle
	Each	Aerobic Blood Culture Bottle
	Each	Pediatric Blood Culture Bottle
	Each	Blue Top Single Swab
	Each	Blue Top Double Swab
	Each	Red Top Single Swab
	Each	Red Top Double Swab
	Each	ESWAB Swabs
	Each	M4 Viral Swab
	Each	Stool Culture Media
	Each	O&P Kit
	Each	Castille Wipes
	Each	Urine Culture Kit
	Pack	Urine Culture Port-A-Cul 5ML Vial
Other		
QTY	UNITS	DESCRIPTION
	Each	Transfer Tubes White or Amber (circle one or both)
	Each	Transfer Tube Caps White, Purple, Amber (circle one or more)
	Each	Transfer Pipettes
	Each	Fit Occult Blood (FOBT)
	Each	Cytology - ThinPrep Vial
	Each	Xpert CT/NG DNA Kit Urine or Swab (circle one or both)

3 Convenient Ways to Order

1. Complete order form and fax to: (845) 437-3158
2. Complete order form jchiumen@health-quest.org
3. Call Customer Service 845-554-1115

All orders will be processed within 48-hours.

Supplies for Vassar Diagnostics Specimens Only